

Temp-Taker[®]

Software/Firmware Release Notes 2020-08-28

This page contains release notes for the Temp-Taker[®] software application and handheld device firmware. The notes are presented in descending date order, with the latest listed first.

v5.0.55 – August 28, 2020

This release includes updated device firmware code.

Changes & Enhancements

- Created new battery-conservation feature that extends battery life 3-to-5 times (depending on device usage). The v5.0.55 software-update automatically enables the feature. The feature may be disabled/enabled via the Settings section of the System → Device tab. If enabled, the Bluetooth radio inside the device is automatically powered off when the device powers down. When the device is turned on, the Bluetooth radio is automatically turned back on. When this battery-conservation feature is enabled, the Temp-Taker[®] app cannot turn on the device remotely because there is no Bluetooth connection.
- Enhanced Bluetooth error-state detection; displays message when there is 'No Active Bluetooth Radio Available', an 'Incompatible Bluetooth Radio' is detected, or an 'Incompatible Windows Version' is running.

Fixes

- Fixed Bluetooth connection issue (that is caused by an error in the Windows Bluetooth driver) by implementing a workaround in the Temp-Taker[®] software code. Eliminated the need to disable the Windows Bluetooth power-saving property, which becomes inaccessible in the Windows feature update, version 2004.
- Fixed vulnerability that prevented the device from engaging auto-off when the host PC is put to Sleep under certain conditions.
- Fixed cause of application error that can occur on startup if the Bluetooth radio in Windows settings is turned off.

v5.0.51 – June 25, 2020

This release includes updated device firmware code.

Changes & Enhancements

- Added reminder alarm weekday settings feature which allows the user to specify the weekdays that Program reminder alarms are enabled.
- Made application compatible with Windows Lockdown (Multi-App Kiosk) mode.
- Added protection against inadvertent column sorting in report view. Before sorting is permitted, the corresponding Group Related Events box must be unchecked in Reports → [gear icon] tab.
- PIN resets, performed via the software app, are now recorded in the Security Log section of the System → Security tab.

Fixes

- Adjusted the 'User' column width when printing to take into account a wider custom 'User' title.

v5.0.41 – March 18, 2020

This release includes updated device firmware code.

Changes & Enhancements


- Increased data-upload throughput.
- Added Hourly, 5 min, and 1 min interval options to the Auto-Import feature.
- Added new Auto Flip feature. When enabled, the display automatically flips upside-down according to the unit's physical orientation.
- Added a label to the main window footer containing the count of units in need of new settings after an auto-import.
- Modified auto-generated filenames to minimize the possibility of filename conflicts. Also, all space characters are now replaced with underscore characters.
- Software now automatically displays the Bluetooth section of the System → Device tab on login, if there are no paired devices.
- Enhanced database archiving feature. Backups are now generated as needed throughout the day.
- Improved user calibration process on the handheld device.
- Handheld device now displays 'Pair Bluetooth' option on welcome screen when the device is not paired.

Fixes

- Fixed error in device reminder alarms for Auto-Import and Auto-Report features.
- Fixed error wherein some Program Group and Program fields in Cooking/Cooling data records were empty.
- Addressed issue where the email account Login Password would not bind after saving changes.

v5.0.29 – December 3, 2019

This is the first Temp-Taker® 5 release (post Beta).

 **NOTE:** The Temp-Taker® 5 software can be installed alongside Temp-Taker® 3 or 4. If so done, the installer will copy the data and settings from the previous software version to Temp-Taker® 5 at the time of installation. Although the old software will continue to function, only one instance of Temp-Taker® can run at a time on a given PC. Also, because of the addition of Bluetooth hardware, the Temp-Taker® 5 software application is not compatible with Temp-Taker® 3 & 4 handheld devices. Please contact ITD Food Safety for upgrade options.

New Features, Enhancements, and Changes (over Temp-Taker® 4)

- All communication between handheld devices and the PC Temp-Taker® application is performed via Bluetooth.
- While a Bluetooth connection is maintained, all device generated data is automatically uploaded from the handheld unit to the PC Temp-Taker® application.
- Multiple data report tabs, with their own dedicated scope, can now be created and customized.
- Multiple auto-generated text and email notifications can be configured to report insufficient data-logging activity.
- Auto-Generated Reports (FKA Auto Report) was enhanced to allow *multiple* configurations.
- Application text/button sizes can be scaled to better suit tablet screen sizes.

- The application layout and text were revised to make software usage more intuitive.
- The *error.log* has been replaced with a new *system.log*. When the file grows beyond a certain size it is backed up to a file containing the timestamp in the name and an empty *system.log* file is started.
- Backing up software settings now produces files with *.dt5* file extensions (versus *.dat* files in Temp-Taker® 3 & 4). While *.dat* files can be imported into Temp-Taker® 5, *.dt5* files cannot be imported into older versions of the software.
- Added 2 new View Options for ISO (YYYY-MM-DD) and European (DD-MM-YY or DD/MM/YYYY) date formats
- Added a second buzzer to the handheld device to increase loudness of reminder alarms.
- Added a second LED to the handheld device to improve visual perception of reminder alarms.

Features Removed (from Temp-Taker® 4)

- USB communication is no longer used.

v4.0.15 – August 31, 2018

This release did not include any changes to device firmware code.

Changes & Enhancements


- The “Device Log-In Screen Title” field in the General Settings → Users tab, now controls the User column header string in all data views.
- Set left and top coordinates of print preview window to zero to account for multi-monitor configurations.
- Increased upper limit of the alarm reminder interval in the Process Options → Cooling tab to 540 minutes.

Fixes

- Fixed error wherein deleting a referenced Conditional Question would sometimes cause an error.
- Corrected the rendering of the cells underneath the “Firmware Version” column in the General Settings → Device History tab.
- Fixed error wherein the `C:\ProgramData\Temp-Taker 4\Archived_Data` folder was sometimes not created before the software attempted to access it

v4.0.10 – June 14, 2017

This is the first Temp-Taker® 4 release. It includes corresponding device firmware code used to upgrade compatible Temp-Taker® 3 devices.

 **NOTE:** Temp-Taker® 3 software does not automatically upgrade to Temp-Taker® 4. If you’re running Temp-Taker® 3, please contact ITD Food Safety for information on how to upgrade to Temp-Taker® 4.

New Features and Enhancements (over Temp-Taker® 3)

- Auto Temp feature detects when probe has been inserted into a hot or cold product and automatically initiates a temperature measurement without any user keypad input. Data results are logged into memory (when appropriate). Device also detects when probe has been extracted from the product and, when appropriate, automatically advances to the next item on the list.
- Enhanced predictive temperature measurement algorithm to achieve greater accuracy and improved speed. The new algorithm is dependent on the response time of the probe. Therefore, the handheld device can now be configured (under Main Menu → Settings) either for the standard Fast-Response probe or the Medium (slower-response) probe.

- Cooking and Cooling modes on the handheld device were simplified and designed to work with the Auto Temp feature. The pre-existing tracking features in these modes were retained and can be optionally used on the fly when desired.
- The following Checklist features and enhancements were added:
 - Simplified user interface requiring less button presses; Capability to enable the previous (Temp-Taker® 3) user interface
 - New *programmable question branching* feature: Allows programming the handheld device to ask a new question based on the answer to the current question – typically used to program custom corrective options when answer is undesired; Compatible with all available question types except Custom Multi Select and Date type questions
 - New *additional comment* feature: Allows the addition of comments to questions; When an additional comment is added to a question, “MORE ►” is appended to the question text when displayed on the handheld device; The additional comment is displayed when the user presses the RIGHT keypad button (►)
 - New *audio/visual alert* feature: Questions can be programmed to make the handheld device activate an audio/visual alert – typically used when branching to a ‘question’ acting as a corrective action directive due to an undesired answer
 - New capability to specify that a question’s answer is not logged into memory – typically used for ‘questions’ acting as general instructions to the user
 - New *temperature type* questions: Create questions that are answered by taking a temperature measurement (fully compatible with Auto Temp feature); Conditional questions can be added when temperature measurement is above, below, and/or within the programmable target range
 - Custom type questions include the following new features and enhancements:
 - Single Desired (available in Temp-Taker® 3): Only one answer is selected and the desired answer must be specified at the time the question is created; Undesired answers are flagged when data is uploaded and displayed in the software application; Enhanced with new related features – e.g. the capability of adding conditional branching questions to any or all answers
 - Single Select (new feature): The desired answer is not specified at the time the question is created; Allows for questions that may have more than one acceptable answer; When answering question on the handheld device, only one answer can be selected
 - Multi Select (new feature): User may select one or more custom answers to the question (up to 10)
 - Option to place custom units on either right or left side of the value in numeric type answers – for example: 100^{PPM} versus Bag #9
 - Increased the number of questions that can be downloaded to a device from 250 to 2,000.
 - New Question and Answer data columns available on the software application capable of displaying the full question and multiple-selected answers
- New ‘Verified By’ feature allows supervisors to approve (i.e. verify) custom data reports
- Program custom corrective options to force user to perform an immediate data upload; this enhancement is typically used with a custom data-reporting service to alert management of a critical condition.

Features Removed (from Temp-Taker® 3)

- The View Save Data mode and 21-day activity counter memory on the device were removed.
- Process symbols were simplified.

v3.7.0 – May 21, 2016

This release includes updated device firmware code available only for compatible devices.

Changes & Enhancements

- Question Checklist feature now supports questions with date or numeric answers.
- Custom single-answer Questions can now download to the device.
- The Auto Import tab now allows the saving of an invalid network path after a warning message is displayed. This allows Auto Import customization to be performed on a computer outside the local network.
- Includes updated USB drivers.

Fixes

- Fixed error wherein the last Auto-Report date was not getting properly initialized when utilizing Auto Import to perform the *initial* Auto-Report configuration. This error would, under certain conditions, result in the generation of a long series of empty reports.

v3.5.89 – April 4, 2016

This release did not include any changes to device firmware code.

Changes & Enhancements

- The SFTP auto-generated CSV filename now captures the 'Default Subtitle' field entered in the General Settings → Printing tab when the file *filenamePrefix.txt* (located in the C:\...\Temp-Taker 3\AutoGenCsv folder) contains the character sequence "<Subtitle>".
- When Station Mode is enabled, the Temp-Taker application now launches without displaying the splash screen and controls that normally close the application, now minimize it instead.

v3.5.88 – March 14, 2016

This release includes updated device firmware code available only for compatible devices.

Changes & Enhancements

- The functionality of the Temp-Taker® 3 – SFTP Edition was ported into the General Edition of the software application. SFTP Mode can now be enabled in the General Settings → Auto Report tab. This setting is now saved in the *.dat* file generated in the General Settings → File tab when backing up Temp-Taker settings.
- The SFTP auto-generated CSV filename now captures the first line of non-white-space characters encountered in the file *filenamePrefix.txt* (located in the C:\...\Temp-Taker 3\AutoGenCsv folder). Organizations can use this field (of the CSV file name) as a location code to determine where the file originated.
- The TILT feature received the following enhancements available only on compatible devices:
 - After the TILT Period of a food item has been started, the user may log additional temperature measurements (by pressing the CENTER keypad button in the TILT Processing Screen) while the item's TILT Period is active.
 - When the UP or DOWN keypad buttons are pressed to mark a TILT Item as 'Consumed' or 'Discarded', the device now requests confirmation before the action is taken.
 - Periodically recurring alarms (with jump-to-item option) can be configured (via the software application) to remind users to take additional temperature measurements on food items with active TILT Periods. These reminder alarms can be configured individually for each TILT Group or Program.

- New Station Mode available that can be used in conjunction with ITD Food-Safety Stations. In this mode the application keeps the parent window maximized.
- The software application now allows the saving of Auto-Import and Auto-Export settings even when there are network path errors currently present. This allows for the customizing of Temp-Taker settings on one computer for the purpose of being applied to other computers.
- When the application is maximized, the parent window no longer covers the Windows taskbar.

Fixes

- None.

v3.5.85 – February 4, 2016

This release did not include any changes to device firmware code.

Fixes

- Fixed error wherein the automatic purging of Calendar Programs older than 12 months caused the unintentional purging of some General Program members.

v3.5.83 – January 22, 2016

This release includes updated device firmware code.

Changes & Enhancements

- The software Auto-Update window now includes the version number of the available update and a link to a webpage with the release notes.
- The “Do not overwrite Auto-Import settings” checkbox in the General Settings → File and Auto Import tabs were changed to “Do not overwrite Auto-Import path”. Now, when this box is checked, only the Network Folder and WebDAV URL paths (with accompanying credentials) are protected from being overwritten by the imported settings file. Similarly, the “Do not overwrite Auto-Report settings” checkbox in the General Settings → File and Auto Import tabs were changed to “Do not overwrite Auto-Report path and emails”. Now, when this box is checked, only the Network Folder and WebDAV URL paths (with accompanying credentials) and Email-To Address(es) and Account Settings are protected from being overwritten by the imported settings file. These changes provide headquarters greater control over the settings of their remote locations while still protecting settings that should be unique to each location.
- Removed the Consumed/Discarded option from the device’s TILT Item processing screen when the TILT Period for the item expires.

Fixes

- Fixed firmware error that under rare circumstances caused TILT data to be created with a duplicate event ID resulting in a data upload error.
- Fixed error wherein the Calendar Programs interface in the software application would sometimes not recognize the need to download new programs to the device when moving from one day to the next.
- Fixed error that would not display the notification that Settings had been modified since the last import whenever changes to Report Templates were saved.

v3.5.81 – December 8, 2015

This release did not include any changes to device firmware code.

Changes & Enhancements

- Updated USB drivers to make fully compatible with Windows 10.

v3.5.80 – November 1, 2015

This release did not include any changes to device firmware code.

Changes & Enhancements

- Removed Auto-Import User Interface controls from the File tab and placed them in a new tab called 'Auto Import'.
- The settings import feature was enhanced with options to exclude the settings related to Auto Report and Auto Import. Two new checkboxes 'Do not overwrite Auto-Import settings' and 'Do not overwrite Auto-Report settings' were created. These new controls were implemented for both manual and auto importing of *.dat* files.
- Added checkbox in the General Settings → Security Options tab corresponding to Supervisor access to the Auto Import tab.
- In the General Settings → File tab, the checkbox 'Include the Users list when importing' was renamed to 'Do not overwrite Users list'. Also, the logic corresponding to the checkbox function was inverted.
- Added new checkbox 'Do not overwrite Printing->Default Subtitle' to the File tab (as well as the Auto Import tab) to replace the previous method of controlling this function by clearing the 'Default Subtitle' field in the saved *.dat* file. The Subway Edition renames this checkbox to 'Do not overwrite Printing->Default Subtitle and Franchise #' and includes protection of the 'Franchise #' by checking the box.
- Modified the application's splash screen.

Fixes

- Fixed error in data view grid where filtering on the Programs column would not display values for Calendar Program names.

v3.5.77 – August 20, 2015

This release did not include any changes to device firmware code.

Changes & Enhancements

- Augmented code related to error handling and recovery when dealing with remembered network connections for the Auto-Report and Auto-Import features.
- The file name that is assigned to auto-generated reports now includes an underscore character when the 'Default Subtitle' field is empty in the General Settings → Printing tab. An underscore character is now used instead of the space character that appears *between* the fields used to construct file names for auto-generated reports.

Fixes

- Fixed error wherein a missed auto-report message window was sometimes being displayed erroneously when the Temp-Taker® application was started.
- Fixed error on Data → Report View tab wherein visibility of the Program Group column was dependent on the visibility of the Program column.
- Clicking the 'Undo Changes' button after performing a column sorting operation, now refreshes the grid with original row order.

v3.5.76 – June 9, 2015

This release includes updated device firmware code.

Changes & Enhancements

- The number of characters that can be entered into the 'Mail-To Address(es)' field (in General Settings → Auto-Report tab) was increased from 248 to 750. This increased the number of email addresses that can be entered from 7 to 20 (where each address has 35 characters).
- Before purging old data the software app now creates an XML file of the data (instead of a CSV file). This helps a customer find archived records more quickly when importing the XML file into Excel. The Temp-Taker app still displays data for two full calendar years plus the current year in course.

Fixes

- Fixed error preventing software auto-update from working on some machines.
- Fixed firmware error that occasionally caused the Corrective Option name for Out-of-Range events not to upload from the device to the PC software app.
- Fixed error in Auto Report tab where appending a final comma or semi-colon at the end of a series of email addresses would cause a validation error.
- Removed extraneous underscore characters in Subject line and file name when generating email from print preview window.

This release includes updated device firmware code.

New Features & Enhancements

- New 'Report Templates' feature allow the saving of up to 20 customized report templates based on either Report View or Supervisor View tabs. Reports are customized by hiding/unhiding columns, re-arranging the column order, and applying data filters.
- New 'Auto-Report' feature replaced the Auto-Email feature in v3.5.5. It provides the following options:
 - Generate PDF reports immediately after every data upload from a handheld device.
 - Generate PDF reports every day for the current or previous day's data; time of day can be set by user.
 - Generate PDF reports every fiscal week for the current or previous week's data; week day and time can be set by user.
 - Generate PDF reports every month for the current or previous month's data; day-of-month and time can be set by user.
- Auto-generated PDF reports (described earlier) can be based on Report View, Supervisor View, or *multiple* Report Templates.
- Auto-generated PDF reports may be set to *auto-email* (as it currently did in v3.5.5) **and/or** *auto-save* to a network folder that the current user has permissions to browse and write to. This also applies to cloud-storage folders (e.g. Google Drive, Dropbox, and OneDrive) installed on a local machine or in a network location. An optional WebDAV URL may be used instead of a network folder.



NOTE: WebDAV is a service that allows files to be shared over the Internet without installing additional software. Your cloud storage provider may provide access to a WebDAV file server.

- New Auto-Update Settings feature allows organizations with multiple locations to **remotely** update the Temp-Taker® settings in all their locations. When this feature is enabled, the Temp-Taker® software app periodically looks for a *.dat* settings file in a predetermined network folder (as configured in the app). If one is found, it compares it with the previously imported file and if it is determined to be different, the new settings file is automatically imported. An optional WebDAV URL may be used instead of a network folder.
- New Custom Answers feature for checklist Questions. Custom answers may now be entered for each Question in Take Temps → Questions tab. Up to 10 different answers may be entered for each Question.



NOTE: The initially selected answer on the device is now always the first answer in the list; this also applies to the 'Yes/No' Question Type.

- New Custom Installer feature allows the inclusion of custom settings (prepared by ITD Food Safety upon customer request) into the software installation process. This allows for custom settings to be loaded automatically upon software installation.
- The following icons on the device now auto-hide when not needed: Log Out, TILT Items, Corrective Sessions, Cooling Items, and Cooking Items. For example, if there are no items in Corrective Session, then the icon is hidden. This makes the interface simpler and easier to navigate. The following changes were implemented in conjunction with this enhancement:
 - The Log In icon is no longer used on the device. Log In is automatically prompted when entering a mode that requires it, e.g. Take Temps. The Log Out icon – if enabled via the software – is only displayed when a user is currently logged in.
 - The positions of the following icons were changed on the device: Spot-Check Temps, Cooling Items, and Cooking Items. The Spot-Check Temps icon was positioned to the left of the Take Temps icon. The Cooling and Cooking icons were positioned to the right of the Take Temps icon.

- The 'Add/Remove Device Features' section of General Settings → Device Settings tab of the software was revised to reflect the changes described above. The following icons were removed from the tab: TILT Items, Corrective Sessions, Take Temps, Cooling Items, and Cooking Items. Since these can no longer be enabled/disabled by the user, they were replaced with "[...]". Tool tips were added to the icons.
- The counter displayed at the top left corner of the TILT Items, Corrective Sessions, Cooling Items, and Cooking Items icons now blink.
- Added a 1-stage hot-food Cooling Process option to the already existing 2-stage hot-food and room-temp food cooling.
- Added option to automatically end the Cooking Process on first In-Range or overcooked temp.
- When the user presses the CENTER key to take the temp of a Holding item, the device now checks to see if the item is currently in Corrective Session. If it is, the user is advised and given an option to jump to the Corrective Session. The same has been implemented for TILT.



NOTE: V3.5.5 already implemented similar functionality for Cooking and Cooling.

- The following enhancements and changes in the Timers mode have been implemented:
 - Custom timer names can now be assigned to each of the 3 Timers via the PC software in General Settings → Device Settings. A maximum of 11 characters can be entered for each name. The name is displayed in the screen title of Timers mode and when a Timer countdown alarm activates.
 - The screen title was revised to make more apparent which Timer is currently displayed.
 - New Auto-Reset feature was added. When Auto-Reset is 'On' the Timer is automatically reset upon counting down to zero and a new countdown cycle restarts.
 - Holding down the CENTER key now displays an Options screen for the corresponding Timer. Options provided are Reset, Clear, and Set Auto-Reset. Selecting 'Reset' restores the latest Timer setting entered before a countdown was started. Selecting 'Clear' zeros-out the Timer. Auto-Reset can be turned on/off by selecting the 'Set Auto-Reset' option.
 - Timer status indicators were added to the bottom of the screen in Spot-Check Temps mode, Timers mode, Clock mode, and the Main Menu.
 - A different audio/visual indication pattern was assigned for each of the 3 Timers. The alarm duration time for each Timer was extended to 1 minute.
 - The countdown settings for Timer 1, 2, and 3 are now being saved in non-volatile memory; as a result, these settings are now preserved when power is removed.
- Added the current time and battery level icon to the bottom of the screen in Spot-Check Temps mode. Also added audio feedback to Spot-Check Temps statistics mode. User now receives both visual and audio indications whenever a new MAX or MIN is detected.
- New 'Quick Help' device feature now periodically displays a Quick Help message at the bottom of the screen to help new users become acquainted with how the interface works. This feature may be disabled at any time in Settings → 'Quick Help Messages' directly on the device.
- TILT Reminder Alarms now activate for every item in the Program or TILT Group. Also, the TILT Period maximum was extended from 17 hours to 1 week.
- Device now prevents the Auto-Reduce Brightness battery-conservation feature (set via the software) from reducing the brightness of non-transflective LCDs to below 5% (unless the LCD brightness level was already set below 5% from within the device's settings). This allows maximizing battery conservation by setting the Auto-Reduce Brightness to 100% *even when* connecting devices with different LCD technologies to the same PC.
- Integrated new Auto-Report features to existing Subway-Report feature in Subway Edition of the software.

Fixes

- Fixed the following error in Timers mode: If user stopped and resumed a countdown timer, the reset setting of the timer would be changed to the time remaining when user resumed the countdown. Now resetting the timer will always restore the original time that was set before timer countdown began.
- Fixed error that prevented device's Auto-Off feature from activating when in the 'Confirm Cooking End' and Cooling 'Options' screens.
- Fixed error wherein importing a settings file that has Programs of type Multi with defined Checklists into a pre-version 3.5 app forces user to reassign all Questions to that Checklist.
- Fixed error that sometimes resulted in missing the Program Group name during a data upload.
- Fixed error wherein replicating a Program or Calendar Program of type TILT or Cooling Group performs a Celsius to Fahrenheit conversion on all temperature fields.
- Fixed error wherein importing a 'Read-only' *.dat* settings file caused future imports to fail.
- Fixed error wherein changes to the Station Functionality section in General Settings → Security Options tab did not persist after restarting the application.

v3.5.5 – December 22, 2014

This release includes updated device firmware code.

New Features & Enhancements

- New software update feature alerts user when a future software version becomes available for download via the internet; software can be easily updated while preserving current settings and data.
- New auto-email feature allows automatically sending a PDF report of the Supervisor View to specified email address(es) after every device-to-PC data upload; when used in conjunction with a 3rd party service (e.g. EmailItIn.com) it can be used to automatically upload reports to a variety of cloud storage accounts.
- New Cooling feature automates tracking of the cooling process for hot & room-temperature food items following HACCP guidelines; provides reminder alarms and corrective option features.
- New Cooking feature automates temperature tracking of the cooking process.
- New Advanced Program Options dialog (accessible in the Program Options section of the Programs tab) includes options to enable the creation of Holding Groups (formerly called Temp Groups), Cooking Groups, Cooling Groups, TILT Groups, and Question Groups (aka Checklists); checking any of these options will cause a corresponding tab to become visible where the group(s) can be created. The dialog also contains an option to 'Download Process Symbols' which determines whether a special symbol character identifying the process (i.e. Holding, Cooking, Cooling, TILT, and Question) is prepended to all names where appropriate.
- The 'Type' column of the Programs grid has been renamed 'Process' and includes new options for TILT, Question, Cooking, and Cooling; choosing any of these options allows adding existing items *directly* into the program for that process. In addition, a 'Multi' option (formerly called Groups) is conditionally offered if Multi-Process Programs are enabled via the Advanced Program Options dialog.
- A new section called Process Options (accessible via the vertical toolbar) contains dedicated tabs for managing the *global* settings of all process types, including the options and zones related to the Holding and TILT processes that were previously found in Temps Setup.
- The optional 'All Items' program can be configured to include either all existing holding zones or all assigned holding zones with each item downloaded to device.
- The interface for creating Checklists (aka Question Groups) is now on a separate tab from the interface for creating Questions, and is consistent with the interfaces for creating Holding Groups and TILT groups.

- Clicking the footer area of the application when the application is maximized on the screen will restore the application to its previous size.
- Improved touch support for reordering rows on editable data grids.
- Improved touch support for displaying right-click context menu to show/hide columns on Temps Report.
- Calibri font now included with installer package for consistent report rendering on XP systems that do not have Microsoft Office installed.
- User is now alerted on device when entering an incomplete Program (i.e. a Program that has been partially checked off).
- Simplified the user directives that appear at the bottom of the device screen.
- Added visual feedback at application launch to indicate that app is searching for connected devices.
- Eliminated Item Categories.
- Eliminated separate view settings for Report View; global view settings now apply.

Fixes

- Resolved error where Zones grid would not allow the re-ordering of rows after the insertion of a new row unless changes were saved and the tab content reloaded.
- Resolved error where wrong row was replicated when choosing 'Replicate Row' from the context menu of a grid while a drop-down was open on another row.
- The Calendar Programs interface is now refreshed if visible while cycling from one day to the next so as to accurately display the current load status.
- The dropdown list of Program Groups in the Programs grid (when Program Groups are enabled) is correctly refreshed after importing a new file.
- Attempting to enable All Items program with a blank name will no longer prevent other valid global settings from being saved.
- Resolved error where deleting a Holding Group (formerly called Temp Group) or TILT Group would occasionally result in removal of an item from a Calendar Program.
- Clicking Undo in the Programs tab will now correctly restore all items within Programs that were deleted prior to the Undo operation.
- Corrected error where the setting to require supervisor PIN for corrected items was not being sent correctly to the device in certain cases.
- The auto-compose Email feature of the Print Preview window has been re-implemented to avoid occasional errors related to multiple Outlook processes.
- Resolved device firmware error where the automatic Program Group selection made during execution of jump options was inadvertently changed by code that automatically selected the Program Group that represented a calendar day with today's date.

v3.1.22 – September 28, 2014

This release includes updated device firmware code.

New Features & Enhancements

- The 'Auto-Load Programs' dropdown in Temps Setup → Calendar Programs tab was enhanced with two new choices: 'Today + 6 days' and 'Today + 30 days'. Whenever one of these choices is selected, connecting the device will download the setup for the current day plus the following 6 or 30 days respectively.
- Added the battery level icon to the bottom of the device screen in Timers and Clock modes.
- Added the current time reading to the bottom of the device screen in Timers mode.

- In Timers mode, the device's UP/DOWN keys now increment/decrement the time in 1-minute steps (previously 15-second steps); the fast step size was increased from 1 minute to 5 minutes.
- Whenever Programs corresponding to more than one calendar day have been downloaded to the device (using the Calendar Programs interface), the handheld device will now automatically make the 'Program Group' corresponding to today's date the default selection every time that Take Temps mode is entered. For example, if all the Programs for the month of August have been downloaded to the device and today is August 8, the handheld device will automatically make the 'Program Group' corresponding to August 8 the default selection whenever Take Temps mode is entered.
- *AutoExport.exe* and *AutoImport.exe* command-line tools now return an exit code to facilitate error detection in a scripted environment.
- Eliminated unnecessary download of entire setup to connected devices when only the Users list has changed.
- Device's Auto Off feature is now prevented from activating if in Timers mode and at least one of the Timers is currently counting down.
- Switching from Temps Data to Temps Setup or Settings when the scope is set to Day and then back to Temps Data will no longer require the user to re-navigate to the desired day, as the previously selected date is remembered.
- Implemented enhancements to custom Subway reporting feature.
- Added the Station Functionality section in the Settings → Security Options tab; this section allows defining the PC on which the software is running as a 'Read/Write Station' (i.e. uploads device data to PC *and* downloads new setup to device), 'Read-Only Station' (uploads device data to PC *only*), or a 'Write-Only Station' (downloads new setup to device *only*).
- Added code to protect data uploaded to a PC that had its system date/time tampered with; data will be output to a CSV file.

Fixes

- Software no longer shuts down when attempting to print data reports displaying the Program Group column if the data set contains records created in an older version of the program where program groups did not exist.
- Software no longer shuts down when encountering a printer utilizing a legacy driver incompatible with Microsoft's .NET Framework.
- Firmware updates are now prevented from occurring on read-only stations (i.e. when the Station Functionality is set to 'Upload temperature data only').
- When the device automatically powers off while on the Main Menu with the Log Out icon selected, the Log In icon is now correctly displayed when at a later time the device is turned on.
- The device's manual reset option – for resetting all checked off items in a Program – is now correctly being provided when the RIGHT key is pressed for Programs containing members of type Group.
- Print Preview and PDF reports now correctly adjust column widths when columns are reordered or when Report View uses different Date/Time formats from Supervisor View; also, previously-missing 'Group #' and 'Seq #' column values are now displayed.
- Resolved error where downloading a large number of question Checklists (>250) would sometimes cause an error message to display on the handheld device when attempting to access certain Programs.
- Resolved problem of device displaying a rounded version of a temperature reading while using the more exact actual measurement behind the scenes when evaluating whether the measurement was In Range or Out of Range.
- Resolved error involving misaligned column widths in Corrective Options grid under certain conditions.
- Clicking the box to disable categories while the Category cell of the Items grid is being edited no longer causes the application to shut down.

- Create PIN dialog is now properly shown after importing a file involving the use of security PINs.
- Name of logged in user in main window is now properly refreshed when user name is modified in Settings.
- Corrected error involving misaligned column widths in Corrective Options grid under certain conditions.
- Clicking the box to disable categories while the Category cell of the Items grid is being edited no longer causes the application to shut down.

v3.1.12 – May 1, 2014

This release includes updated device firmware code.

New Features & Enhancements

- New custom Subway reporting feature.
- New column added to Items grid which allows selecting between Probe and Keypad input methods. When Keypad input is selected the device will request the temperature measurement to be inputted via the keypad; this feature is useful when 'taking the temp' of a cooler or freezer equipped with an accurate external temperature display.
- New checkbox added to Corrective Options tab that enables requiring a supervisor PIN when recovering items in corrective session. When this box is checked, the device will request a supervisor PIN be entered when an In-Range temp is taken on an item in Corrective Session. This new capability was implemented to make Temp-Taker® 3 fully compliant with Subway operation guidelines.
- Improved device firmware execution speed by optimizing functions used to find objects.
- Improved USB detection, including the detection of connected devices on application launch.
- Improved error detection and UI responsiveness when non TT3-devices are connected via USB.
- Additional time was added before performing a device firmware update when drivers are loaded for the first time; this change avoids driver interference with firmware update process.
- *AutoImport.exe* command-line tool now includes help switch and recognizes import file when extension (.dat) is omitted.

Fixes

- Software now properly decodes the degree character (°) when uploaded back from the device.
- Cancelling out of the Open File dialog used to import a setup file no longer blocks subsequent communication sessions when device is connected.
- Importing a setup file that requires creating an admin PIN no longer causes the backup file to appear as not current.
- Viewing older data collected from software versions prior to 3.1.0 no longer causes the application to shut down when the Program Groups column is visible.

v3.1.0 – March 3, 2014

This release includes updated device firmware code.

New Features & Enhancements

- New optionally-enabled Program Groups allow aggregating Programs in collections that can be selectively downloaded to specific devices identified by the serial number.
- Severed interdependency between General Programs and Calendar Programs with regard to program load counts. Calendar interface now allows loading a maximum of 15 programs per day.

- As with general programs, calendar programs can now be created specifying a member type of 'Groups' or 'Items', the latter of which allows adding individual Items assigned to any zone directly in the Program.
- New *AutoImport.exe* and *AutoExport.exe* command line tools allow creating scripts to automatically import setup changes or export data reports with specific parameters without requiring direct user interaction.
- Can now create reminder alarms for calendar programs.
- New Auto-Load feature for Calendar Programs allows enabling the software to automatically load all programs in a given scope (day, week, or month) when device is connected. Calendar programs' load status is automatically refreshed when cycling from one day to the next.
- Can now copy/paste entire day or entire week of calendar programs via right-click context menu.
- Added special formatting to easily identify empty calendar programs.
- Users tab now allows customizing the Log-In screen title on the device.
- Improved location placement for message dialog windows.
- The device now skips over single-element listings to save button presses whenever possible.

Fixes

- Resolved error on device firmware that caused some reminder alarms to be missed.
- Resolved error on device firmware that on rare occasions caused the LCD backlight to stay on after 'powering down' the device.
- Help dialog of PC application now ensures the software version number is always visible to users across systems.
- Corrected the calculation of column widths when printing, previewing, or exporting Temps Data reports.
- Column sorting indicators are now reset after Undo or Save operation on all grids.
- Checklist dialog no longer attempts to delete empty parent program when deleting a checklist.
- USB communication with connected device is now blocked if a setup import operation is in progress.

v3.0 RTM – April 11, 2013

This was the first Temp-Taker® 3 release.